TCL200

User Guide

Digital Alarm Clock Features	3	Clock Functions	12
Cordless Vibration Pad Features	5	Help and Support	16
Important Safety Instructions	6	Troubleshooting Guide	16
Included Equipment.....	7	Maintenance and Warranty.....	17
Setup	7	Technical Details.....	19
LED Display Notification.....	11		

FCC Compliance


The FCC requires that you will be advised of certain requirements involving the use of this telephone.

1. This equipment complies with Part 68 of the FCC rules. A label on the base unit of this equipment contains, among other information, the ACTA Registration Number (US:AAAEQ##TXXXX) and the Ringer Equivalence Number (REN) for this equipment. You must, upon request, provide this information to your telephone company.
2. An applicable certification jacks Universal Service Order Codes (USOC) for the equipment is provided (i.e., RJ11C) in the packaging with each piece of approved terminal equipment.
3. A plug and jack used to connect this equipment to the premises wiring and telephone network must comply with the applicable FCC Part 68 rules and requirements adopted by the ACTA. A compliant telephone cord and modular plug is provided with this product. It is designed to be connected to a compatible modular jack that is also compliant. See installation instructions for details.
4. The REN is useful to determine the quantity of devices you may connect to your telephone line and still have all of those devices ring when your number is called. In most, but not all areas, the sum of the REN of all devices should not exceed five (5.0). To be certain of the number of devices you may connect to your line, as determined by the REN, you should call your local telephone company to determine the maximum REN for your calling area.
NOTICE: If your home has specially wired alarm equipment connected to the telephone line, ensure the installation of this equipment does not disable your alarm equipment. If you have questions about what will disable alarm equipment, consult your telephone company or a qualified installer.
5. If your telephone causes harm to the telephone network, the telephone company may discontinue your service temporarily. If possible, they will notify you in advance. But if advance notice is not practical, you will be notified as soon as possible. You will be advised of your right to file a complaint with the FCC.
6. Your telephone company may make changes to its facilities, equipment, operations or procedures that could affect the proper operation of your equipment. If they do, you will be given advance notice so as to give you an opportunity to maintain uninterrupted service.
7. If you experience trouble with the telephone, please contact the Amplicom Service Center at 1-866-AMPLICOM (267-5426) for repair/warranty information. The telephone company may ask you to disconnect this equipment from the network until the problem have been corrected or you are sure that the equipment is not malfunctioning.
8. This equipment may not be used on coin service provided by the telephone company. Connection to party lines is subject to state tariffs. (Contact your state public utility commission or corporation commission for information.)

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation. Privacy of communications may not be ensured when using this phone. Changes or modifications not expressly approved in writing by Amplicom may void the user's authority to operate this equipment.


Some cordless phones operate at frequencies that may cause interference to nearby TVs and VCRs; to minimize or prevent such interference, the base of the cordless phone should not be placed near or on top of a TV or VCR; and, if interference is experienced, moving the cordless telephone farther away from the TV or VCR will often reduce or eliminate the interference.

TCL200 Digital Alarm Clock Features

Front View


- 1 Snooze button /Voice prompt /Incoming call light /Stop alarm
 - 2 LED: Alarm 1 ON/OFF
 - 3 Alarm ON/OFF
 - 4 LED: Alarm 2 ON/OFF
 - 5 Repeat alarm signal
 - 6 Increase value /Clock voice prompt
 - 7 Decrease value /Select language
 - 8 Set weekday and time
 - 9 Set alarm time and alarm activation period
 - 10 LED: Alarm period, Monday to Friday ON/OFF
 - 11 Display
 - 12 LED: PM (for 12-hour clock)*
- *Only when English is selected

Back View

- 13** Volume control
- 14** Set one of 5 acoustic signals
- 15** Pitch
- 16** Vibration and Acoustics
- 17** Flashing lamp ON/OFF
- 18** Snooze delay
- 19** Dimmer ON/OFF
- 20** Telephone line socket
- 21** Power adapter socket
- 22** Vibration pad socket
- 23** Battery compartment
- 24** Set alarm period

TCL200 Cordless Vibration Pad Features


1 LED: Alarm ON/
OFF/Battery
charge low

2 Vibration pad
ON/OFF

3 Battery
compartment

4 Screws
for battery
compartment

5 Connection
cable socket

Important Safety Instructions

Read this operating manual thoroughly.

SAVE THESE INSTRUCTIONS

Intended use

This alarm clock is designed to provide an acoustic alarm or alert signal at a set time or on receipt of a telephone call. Any other use is considered unintended use. Unauthorized modification or reconstruction is not permitted. **Do not** open the device and attempt any repair work yourself.

Installation location

Prevent excessive exposure to smoke, dust, vibration, chemicals, moisture, heat and direct sunlight.

Heart pacemaker

Always maintain a distance of at least 6 inches from implanted heart pacemakers to prevent any possible interference.


Getting Started — Included Equipment

- 1 Digital Alarm Clock
- 1 Power Adapter
- 1 Vibration Pad
- 1 Connector Cable for Vibration Pad
- 2XAAA, 500mA Ni-MH Rechargeable Battery for Alarm Clock
- 3XAAA, 800mA Ni-MH Rechargeable Battery Pack for Vibration Pad
- Telephone Cord
- Operating Instructions
- Battery Compartment Cover

Setup

Connecting the alarm clock

1. Insert the supplied batteries in the clock's battery compartment.


Important: Pay attention to correct polarity when inserting the batteries.

2. Replace the battery compartment cover on the underside of the alarm clock.
3. Screw the battery compartment cover to the alarm clock housing.

Connect vibration pad and phone cord

Connect the equipment as illustrated below:

Important: To connect the alarm clock to a telephone connection for use as a signal amplifier, you may need an appropriate adapter.

Important: Only use the rechargeable battery pack model GPHCH73N07, 800mAh NiMH supplied by GPI INTERNATIONAL LIMITED and the telephone cord provided.

Set the language

1. Press the snooze button to check the current language. The alarm clock announces the current time and weekday.

2. To change the language, hold down the select language button for 10 seconds. The alarm clock will announce the time and weekday — in the next language.
3. The languages will change in the following order: German, French, Italian, Spanish, and English.

Set the weekday and time

1. Press the Set Weekday and Time button for 5 seconds. The alarm will announce the weekday currently set.
2. Use the Increase/Decrease Value buttons as needed to set the current weekday. The alarm clock will announce the weekday currently set.

Important: If the language is set to English, you can then select the time format. Press the Increase/Decrease Value buttons to select between the 12 and 24-hour display. Press the Set Weekday and Time button to confirm the setting. Skip Step 3 below.

3. Press the Set Weekday and Time button to confirm the setting. The alarm clock will announce the time currently set.
4. Use the Increase/Decrease Value buttons as needed to set the current hour, and then press the Set Weekday and Time button. The alarm clock will announce the minute currently set.

5. Use the Increase/Decrease Value buttons as needed to set the current minute, and then press the Set Weekday and Time button. The alarm clock will announce the complete weekday and time setting.
6. Clock voice prompt on/off setting: Press and hold “+” button for 10 seconds. If the voice prompt is weekday and time, the setting will be on. If the voice prompt is two beeps, the setting is off.

Setting the alarm time and period

Important: There are two alarm period settings available. You can set one of three alarm periods for each setting (Monday to Friday, Monday to Sunday, or Saturday and Sunday). The alarm time you set is repeated only on the days set for the selected alarm period.

1. Use the Alarm Period switch to select an alarm period (either Alarm Period 1 or Alarm Period 2).
2. Press the Set Alarm Time button for 5 seconds. The alarm clock will announce the currently set hour.
3. Press the Increase/Decrease Value button as needed to set the hour for the alarm, and then press the Set Alarm Time button. The alarm clock will announce the currently set minute.
4. Press the Increase/Decrease Value button as needed to set the minute for the alarm, and then press the Set Alarm Time button. The alarm clock will announce the alarm period currently set.

5. Press the Increase/Decrease Value button as needed to set the alarm period, and then press the Set Alarm Time button. The alarm clock will announce the complete alarm setting. The alarm volume/acoustic/pitch settings are available in case of the alarm on/off switch is ON and alarm period switch is setting to 1 or 2 but not period 1+2.

LED Display Notification

Vibration Pad

The following functions are on the vibration pad:

Flashing Green Every 2 Seconds: Vibration pad is on.

Red: Vibration pad is charging.

Flashing Red: Battery charge is low.

Clock Functions

Certain functions require entering test mode before setting.

Turning the alarm on/off

Turn the alarm on/off rotary knob clockwise for On. Counterclockwise for Off.

Setting the volume

Caution: The alarm and telephone ringer begin sounding at a low volume. The volume will then increase to the set maximum until the alarm is turned off or the telephone is answered.

Turn the volume rotary knob to set the maximum volume for the alarm and telephone ringer. Set the volume to a quiet setting to start and increase it slowly to the desired maximum volume.

Caution: Exposure to the alarm at full volume for a long period could cause ear damage.

For test mode

Turn the alarm on/off switch to the on position. Turn all switches on the back of the clock to their far right position (when looking at the back of the clock). Press and hold the Snooze button and plug in the clock. Press the snooze button again to activate Test Mode.

Selecting the acoustic signal


Clock must be in Test Mode before selecting the acoustic signal. Turn the acoustic signal rotary knob to required acoustic signal.

Setting the pitch

Clock must be in Test Mode before setting the pitch. Turn the pitch rotary knob to set the pitch. It may be possible to detect higher or lower pitch tones even at low volumes.

Setting the clock face dimmer

Slide the dimmer switch to activate clock face dimmer. When set to ON, the brightness of the clock face is automatically adapted. When set to OFF, automatic adaptation is switched off.

Selecting the alarm period


Slide the alarm period switch to select between the alarm periods. There are three settings available: (from left to right) alarm period 1, alarm period 2, or alarm periods 1 + 2.

Selecting the alarm mode

Slide the alarm mode switch to set the type of alarm. There are three settings available: (from left to right) vibration only, acoustic signal only or acoustic signal and vibration.

Setting the flashing light

Slide the switch to set the light in the snooze button to flash during an alarm. There are two settings available: (from left to right) On and Off.

Setting the snooze delay


Slide the snooze delay switch to set the snooze delay. There are three settings available: (from left to right) 0 min, 5 min, or 15 min.

Activating the snooze delay

Press the snooze button to stop the alarm for a specified time before it resumes. To activate the snooze function, see "Setting the snooze delay".

Activating the vibration pad

After connecting the vibration pad, use the switch on the vibration pad to turn the vibration function on and off. There are two settings available: (from left to right) On and Off.

Setting the alarm time and alarm period

Press the set alarm time and alarm period button to set the alarm time and the active alarm period.

Increase selection value

Press the increase value button to increase the selected value in the display.

Decrease selection value

Press the decrease value button to decrease the selected value in the display.

Repeat alarm

Press the repeat alarm button while the alarm is sounding to switch off the alarm and set the same alarm time for the following day.

Help and support

If you have problems with your device, please refer to the following information before calling for support. If you continue to have problems, please call the Service Center at 1-866-AMPLICOM (267-5426) for further assistance.

Troubleshooting guide

My clock does not work

1. The battery charge may be too low or no batteries have been inserted.
2. The clock may not be plugged in.

My vibration pad does not work

1. The battery charge may be too low or the battery is defective (additional Battery packs are available from Amplicom).
2. The vibration pad may not be switched on.
3. The vibration pad is too far away from the alarm clock.

The alarm does not go off

The alarm function is not switched on; turn the knob to “On”.

Maintenance and Warranty

Maintenance

- Please clean your equipment's surfaces with a soft, lint-free cloth.
- Never use cleaning agents or solvents.

Warranty

AMPLICOM equipment is produced and tested according to the latest production methods. The implementation of carefully chosen materials and highly developed technologies ensure trouble free functioning and a long service life. The terms of the warranty do not apply where the cause of the equipment malfunction is the fault of the telephone network operator or any interposed private branch extension system. The terms of the warranty do not apply to the rechargeable battery pack or power packs used in the products. The period of warranty is 24 months from the date of purchase.

All deficiencies resulting from material of production faults which occur during the period of warranty will be eliminated free of charge. Rights to claims under the terms of warranty are annulled following intervention by the purchaser or third parties. Damage caused as a result of improper handling or operation, incorrect positioning or storing, improper connection or installation, Acts of God or other external influence are not covered by the terms of warranty.

In the case of complaints, we reserve the right to repair or replace defect parts or provide a replacement device. Replacement parts or devices become our property.

Rights to compensation in the case of damage are excluded where there is no evidence of intent or gross negligence by the manufacturer.

If your equipment shows signs of defect during the period of warranty, please return to the sales outlet in which you purchased the AMPLICOM equipment together with the purchase receipt. All rights to claims under the terms of warranty in accordance with this agreement must be asserted exclusively with regard to your sales outlet.

Two years after the purchase of our products, claims under the terms of warranty can no longer be asserted.

Declaration of conformity

The ETL Listed Mark is an alternative to the CSA and UL marks. ETL Testing Laboratories owned by Intertek Testing Services (ITS) is recognized by OSHA as a Nationally Recognized Testing Laboratory (NRTL). ITS tests products according to nearly 200 safety and performance standards. The ETL Listed Mark is accepted throughout the United States when denoting compliance with nationally recognized standards such as ANSI, IEC, UL, and CSA.

This certification mark indicates that the product has been tested to and has met the minimum requirements of a widely recognized (consensus) U.S. product safety standard, that the manufacturing site has been audited, and that the applicant has agreed to a program of periodic factory follow-up inspections to verify continued conformance.

Technical details

Alarm Clock

Frequency	433 MHz
Power supply	Input: 210-240 VAC, 50/60Hz Output: 12VDC, 670mA
Batteries	2 x AAA 1.2 V, 500 mAh, NiMH
Initial charging time	Up to 16 hrs
Optimum ambient temperature	32°F to 104°F (0°C to 40°C)
Optimum storage temperature	-4°F to 140°F (-20°C to 60°C)
Transmission power	10 dBm
Dimensions L / W / H	3.1 x 5.1 x 5.8 in (80 x 130 x 146 mm)
Weight	12.2 oz (346 g)

Vibration pad

Range	Approx. 328 ft (clear space)
Batteries	3 x AAA 3.6 V, 800 mAh, NiMH
Initial charging time	Up to 16 hrs
Optimum ambient temperature	32°F to 104°F (0°C to 40°C)
Optimum storage temperature	-4°F to 140°F (-20°C to 60°C)
Dimensions L / W / H	4.5 x 4.3 x 1.4 in (115 x 110 x 35 mm)
Weight	6.3 oz (178 g)

THIS PAGE LEFT INTENTIONALLY BLANK

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Distributed by: Amplicom USA LLC
Distribuido por: Amplicom USA LLC