

Geemarc Ampli200+

EASY-DIAL PHONE WITH VOLUME AND TONE ADJUSTMENT



Caution: This telephone product is specifically designed for people with hearing difficulties and is provided with a handset earpiece volume control for individual requirements. Due care must be taken by all users that the handset earpiece volume control is set to the lowest level acceptable by each respective user. Care should therefore be exercised to ensure that any new users are familiar with this requirement.


English

Español Français

CONTENTS

DESCRIPTION	3
SETTING UP	
- Connecting the telephone	4
- Alarm and ring options	4
- Ringer settings	4
- Tone/Pulse switch (T/P)	5
- Time break recall switch (600/100)	5
- Hotel message waiting switch	5
- Wall mounting	6
USING THE PHONE	
- Answering incoming calls	7
- Volume and tone receiving	7
- Making a call	8
- Last number redial	9
- Recall, * & # buttons	9
- Mute button	9
- Storing telephone numbers	10
- Dialling stored numbers	11
REGULATORY COMPLIANCE	12
GUARANTEE	14

DESCRIPTION



(*) : For these optional accessories, consult the retailer (see contact details at the end of this manual)

SETTING UP

1 - CONNECTING THE TELEPHONE

- Connect the handset cord
- Open the battery door by unscrewing the retaining screw. Insert 4 x AAA alkaline batteries. Snap the battery door in place.



In case of power failure, batteries are required for the ring.

- Connect the line cord .
- Connect the PTT plug to the wall socket.
- Connect the 12V DC adaptor cord to the telephone.
- Connect the 12V DC adaptor to your wall socket (**).

(**) Classified “hazardous voltage” according to EN60950 standard.

2 - ALARM AND RINGER OPTIONS

A 3 position switch located on the right side of the phone allows the settings below:

- **0** : Only the ringer is activated.
-  : Ringer and shaker are activated.
-  : Ringer and strobe flash are activated.

3 - RINGER SETTINGS

A switch located on right side allows the ringer to be switched HI,LOW,OFF:

 : Off,  : Low,  : Hi

Your phone comes with the ringer switched to «High» and the melody set to «6».

SETTING UP

	Low	Middle	High
Melody 1	1	2	3
Melody 2	4	5	6 (by default)
Melody 3	7	8	9

If you don't like the ringer melody, you can choose a new one by changing the settings as follow:

- Lift the handset.
- Press **PROG** button.
- Press **#** button.
- Press any button 1 to 9.

Ringer will be activated on next call.

4 - TONE/PULSE SWITCH (T/P)


The Tone/Pulse switch is factory pre-set to Tone position (T). In the US, all telephone exchanges now use Tone dialling. If your phone does not dial out, it is probably being used from an older private switch-board (PBX). In this case, the Tone/Pulse switch must be moved to Pulse.


5 - TIME BREAK RECALL SWITCH (600/100)

For the US, the TBR switch is factory pre-set to 600ms. It can be switched to 100ms for other locations or PBX's.

6 - HOTEL MESSAGE WAITING SWITCH

These are used with the new services provided on hotel network.

A switch  located behind your phone allows to set **ON** the hotel message waiting function.

When the switch is setting **ON** and when you receive a message, the red light  is lit to inform you that you have a message on phone .

SETTING UP

7 - WALL MOUNTING

To wall-mount your telephone, slide out the clip, turn around and slide back into the wall mounting position (located in the base unit, just below the ear piece). This will ensure that the handset stays firmly in position. Make 2 holes in the wall 80mm

apart, insert the wall plugs and the wall mounting screws. The screws should protrude from the wall by 6-7 mm. Place the phone onto the screw-heads and slide down to secure.

If you are using this telephone with a hearing aid, set your hearing aid to the “T” setting.

USING THE PHONE

1 - ANSWERING INCOMING CALL

Without headset:

- When an incoming call is received, the phone will ring, to answer the call lift the handset and speak.
- On completion of the call, carefully replace the handset in the cradle to release the line.

With headset (option):

- When an incoming call is received, and if you wish to answer the call with the headset, press the **SPKR** button and speak.
- On completion of the call, press the **SPKR** button again to release the line.

In handsfree mode:

- When an incoming call is received, the phone will ring, to answer the call press the **SPKR** button.
- On completion of the call, press the **SPKR** button to release the line.



2 - VOLUME AND TONE RECEIVING ADJUSTMENT

You can increase the volume and adjust the receiving tone according to your own hearing requirements. The button **AMPLI** allows extra-receiving and tone functions to be activated or not during the conversation. When the function is activated, the red light of extra-receiving and tone function located near the button **AMPLI** is lit. Adjust the **TONE** and **VOLUME** buttons to remove the level according to your own hearing requirements.

Note:

When the switch **AMPLI** located behind the phone is **ON**, the extra-receiving and tone functions are automatically activated each time you use the phone. When the switch **AMPLI** is switched **OFF**, you will need to activate the function by pressing the button **AMPLI**. Tone and volume adjustment (15dB to 30dB) is only activated when the button is set to **ON**. When this button is set to **OFF** only the volume can be adjusted from 0dB to 15dB

USING THE PHONE



3 - MAKING A CALL

Without headset:

- Lift the handset.
- Await the dial tone and dial required number.
- On completion of the call, carefully replace the handset in the cradle.

With headset (option):

- Press the **SPKR** button.
- Await the dial tone and dial required number.
- On completion of the call, press the **SPKR** button again to clear the line.

USING THE PHONE

In handsfree mode

- Press the **SPKR** button.
- Await the dial tone and dial required number.
- On completion of the call, press the **SPKR** button to release the line.
- If your correspondent doesn't answer your call, press the **SPKR** button again to release the line.

Note:

- To deactivate the handsfree function during a call, lift the handset. To activate the handsfree function during a call, press the **SPKR** button and place the handset on the cradle.
- During a call in handsfree, you can adjust the speaker volume to a desirable level.
Simply use the slide control VOLUME located on the left of the unit.

4 - LAST NUMBER REDIAL

- Lift the handset.
- Await the dial tone and press **REDIAL** button.
- The previously dialled number will be automatically redialled (not in use for memory numbers).

5 - FLASH, * & # BUTTONS

These are used with the new services provided on digital exchanges. For details please contact your network operator.

6 - MUTE BUTTON

If you wish to speak privately to someone else in the room without your caller overhearing, press and hold the button **MUTE**. You will still be able to hear the callers but they cannot hear you during the mute operation. To resume normal conversation, release the **MUTE** button.

USING THE PHONE

7 - STORING TELEPHONE NUMBER

Your phone can store 12 direct memories including 3 priority numbers (21 digits maximum for each memory).

For example, the 3 priority **P1**, **P2**, **P3** memories allow you to store the 3 numbers below:

- Local, Police station,
- Your doctor
- Family

It is recommended to store these numbers during installation

To store a number please follow instructions below:

- Lift the handset.
- Press **PROG** button .
- Press M1, M2, ..., or M9 or one of 3 priority numbers to store your number.
- Dial the number.
- Press button again.

Numbers will be lost after a few hours if you disconnect your phone. Entering a new number in memory will automatically erase the previous number.

USING THE PHONE

8 - DIALLING STORED NUMBER

- Lift the handset and await the dial tone.
- Press appropriate button M1, M2, ..., or M9 or one of the 3 priorities P1, P2, P3 numbers , the stored number will be dialled automatically.

To register your stored phone numbers, use the index card located under the base.

Use this table below to prepare the phone book:

	1	
	2	
	3	
	4	
	5	
	6	
	7	
	8	
	9	
example	Bruce	1 01 64 57 ...

REGULATORY COMPLIANCE

Part 68 of FCC rules information

- This equipment complies with Part 68 of the FCC rules and the requirements adopted by the ACTA. On the bottom of this equipment is a label that contains, among other information, a product identifier in the format US:AAAEQ##TXXXX. If requested, this number must be provided to the telephone company.

- All applicable certification jacks Universal Service Order Codes (USOC) for the equipment is provided (i.e. RJ11C) in the package with each piece of approved terminal equipment.

- A plug and jack used to connect this equipment to the premises wiring and telephone network must comply with the applicable FCC Part 68 rules and requirement adopted by the ACTA. A compliant telephone cord and modular plug is provided with this product. It is designed to be connected to a compatible modular jack that is also compliant. See installation instructions for details.

- The Ringer Equivalence Number (REN) is an indication of the maximum number of terminals allowed to be connected to a telephone interface. The termination on an interface may consist of any combination of devices subject only to the requirement that the sum of the Ringer Equivalence Numbers of all the devices does not exceed five (5). To be certain of the number of devices that may be connected to a line, as determined by the total RENs, contact the local telephone company. For products approved after July 23, 2001, the REN for this product is part of the product identifier that has the format US:AAAEQ##TXXXX. The digits represented by ## are the REN without a decimal point (e.g., 03 is a REN of 0.3).

For earlier products, the REN is separately shown on the label.

- If this equipment causes harm to the telephone network, the telephone company will notify you in advance that temporary discontinuance of service may be required. But if advance notice isn't practical, the telephone company will notify the customer as soon as possible. Also, you will be advised of your rights to file a complaint with the FCC if you believe it is necessary.

- The telephone company may make changes in its facilities, equipment, operation or procedures that could affect the operation of the equipment. If this happens the telephone company will provide advance notice in order for you to make necessary modifications to maintain uninterrupted service.

- If trouble is experienced with this equipment, for repair or warranty information, please contact our customer service staff at: 248-577-5400

If the equipment is causing harm to the telephone network, the telephone company may request that you disconnect the equipment until the problem is resolved.

- Please follow instructions for repairing if any; otherwise do not alter or repair any part of device except as specified.

REGULATORY COMPLIANCE

Opening the equipment or any attempt to perform repairs will void the warranty. For service or repairs, call 248-577-5400

- Connection to party line service is subject to state tariffs. Contact the state public utility commission, public service commission or corporation for information.

If your home has specially wired alarm equipment connected to the telephone line, ensure the installation of this telephone equipment does not disable your alarm equipment. If you have questions about what will disable alarm equipment, consult your telephone company or a qualified installer.

- This equipment is hearing aid compatible.
- This product meets the applicable Industry Canada technical specifications.

GUARANTEE

Sonic Alert warrants the phone against any defect in materials or workmanship for the period of one year from the date of purchase.

Should you experience a problem, contact our customer service department.

Be sure to save your sales receipt as proof of purchase date should you need warranty service.

Within a period of one year from purchase date, Sonic Alert will repair or replace (our discretion) your phone at no cost, if a defect in materials or workmanship is found. If we elect to replace your phone, we may replace it with a new or reconditioned product of the same or similar design.

Repair or replacement will be warranted for a period of 90 days or the original time on the original warranty, whichever is longer.

The warranty does not cover accidents, negligence or breakage to any parts. This includes shipping damage, failure to follow instructions, misuse, fire, floods, use of incompatible accessories, Acts of God or failure in your phone service carrier's line service. The product must not be tampered with or taken apart by anyone who is not an authorized

Sonic Alert representative. Tampering with the phone will void ant written or implied warranties

Sonic Alert shall not be responsible for loss of time, inconvenience, property damage caused by your phone or any other accidental or consequential damages

Warranty service is available only with proof of purchase.

GUARANTEE

Simply send the Ampli200+ (postpaid) and a copy of your sales slip as proof of purchase to:

Sonic Alert Inc., 1050 E. Maple Road, Troy MI 48083

IMPORTANT : YOUR RECEIPT IS PART OF YOUR WARRANTY AND MUST BE RETAINED AND PRODUCED IN THE EVENT OF A WARRANTY CLAIM

SonicAlert

1050 East Maple Rd.

Troy, MI 48083

Ph: 248577 5400

Fax: 248-577-5433

www.sonicalert.com