

BlueHook™ Bluetooth Neck Loop



Wirelessly Connects your Hearing Aid to any Bluetooth™ enabled Mobile Telephone

Product Features

Compatible with Bluetooth™ enabled devices

Bluetooth™ V2.0

Interference-free connection to T-coil hearing aids

Talk time - up to 12 hours

Standby time - up to 250 hours

Charging time - 2-3 hours

Receiving range - 10m

Bluetooth™ frequency 2.4GHz - 2.4835GHz ISM band

3 charging options - USB, Car charger and AC adaptor - all included

One Year Warranty